


How to Add a Button for a Specific Video List

Access Designer from Account's Details Page



Edit Main Menu on Account

The screenshot displays the GrandCare Designer web application. The top navigation bar is dark blue with a hamburger menu icon and the text 'GrandCare Designer'. A left sidebar contains four menu items: 'Accounts', 'Assessments', 'Menus', and 'Logout'. The main content area is titled '327 N Main St Kitchen' and has three tabs: 'MENUS' (selected), 'ASSESSMENTS', and 'SETUP'. Below the title, there are two orange buttons: 'PUSH MENU' and 'ADD'. The 'Menus' section lists several menu items, each with an edit icon (pencil) and a delete icon (trash). A green arrow points to the edit icon for the 'Main Menu' item. A tooltip is visible over the 'Main Menu' edit icon, displaying the text 'Edit this menu on A'.

GrandCare Designer

327 N Main St Kitchen

MENUS ASSESSMENTS SETUP

Menus

PUSH MENU ADD

Main Menu
main - Any Brand

Games
games - Any Brand

Health Menu
health - Any Brand

HEALTHWISE
healthwise - Any Brand

Message Center
message-center - Any Brand

Play Menu
play - Any Brand

Spiritual
spiritual - Any Brand

Edit this menu on A

Add Button to Main Menu

GrandCare Designer

Accounts

Assessments

Menus

Logout

Main Menu

Calendar

Photos

* Name or Description
Main Menu

LAYOUT TYPE: 6 BUTTON

* Menu ID
main

Brand
Any Brand

CANCEL


ADD BUTTON

EXPORT MENU

IMPORT MENU

SAVE MENU

Add Specific Video List



The screenshot shows a mobile application interface with a dark blue header bar containing a grid icon and the text "Add Button". Below the header, the form is divided into several sections:

- Icon:** A small icon of a cooking pot with a flame, followed by the text "Cooking". To the right of this text is a small downward-pointing triangle indicating a dropdown menu. A green arrow points to this triangle.
- * Button Text:** The text "Alex Meal Plan Recipes" is entered into the field.
- * Button Action:** The text "Launch a main screen module" is entered into the field. To the right of this text is a small downward-pointing triangle indicating a dropdown menu. A green arrow points to this triangle.
- * Module Name:** The text "Video List (Specific Video List)" is entered into the field. To the right of this text is a small downward-pointing triangle indicating a dropdown menu. A green arrow points to this triangle.
- * Video List ID:** This field is currently empty.


At the bottom of the form, there is a light orange bar with two buttons: "CANCEL" and "SAVE BUTTON".

View Video List to Grab ID from URL

The screenshot shows a web browser window with the URL `https://my.grandcare.com/content.php?op=viewaccountvideo&listId=5f974dcb-d3e6-4ed8-bda2-7129b5ea5022` highlighted in the address bar. The page title is "Alex Meal Plan Recipes". On the left is a sidebar menu with options: Home, My Profile, Alerts (50), Management, Accounts, Users, Designer, Enterprise, Distributors, Groups, Systems, and Connectors. The main content area displays a modal window titled "ALEX MEAL PLAN RECIPES" with a toolbar containing "Back", "Add Video", "Edit", and "Delete". Below the toolbar is a search bar and a "Refresh" button. A list of three video entries is shown, each with a thumbnail, a title, and a delete icon (X) on the right:

Thumbnail	Title	Action
	On My Own - Living With Autism Let's Make Spaghetti and Meatballs!	X
	On My Own - Living With Autism Coffee	X
	On My Own - Living With Autism Let's Make Avocado Toast!	X

Paste Video List ID and Save



Add Button

Icon
Cooking

* Button Text
Alex Meal Plan Recipes

* Button Action
Launch a main screen module

* Module Name
Video List (Specific Video List)

* Video List ID
5f974dcb-d3e6-4ed8-bda2-7129b5ea5022

CANCEL SAVE BUTTON

The image shows a mobile application interface for adding a button. The dialog box has a blue header with a grid icon and the title 'Add Button'. Below the header, there are four input fields, each with a label and a dropdown arrow. The first field is 'Icon' with a cooking icon and the text 'Cooking'. The second field is '* Button Text' with the text 'Alex Meal Plan Recipes'. The third field is '* Button Action' with the text 'Launch a main screen module'. The fourth field is '* Module Name' with the text 'Video List (Specific Video List)'. Below these fields is a field for '* Video List ID' with the text '5f974dcb-d3e6-4ed8-bda2-7129b5ea5022'. At the bottom of the dialog box are two buttons: 'CANCEL' and 'SAVE BUTTON'. Two green arrows point to the 'Video List ID' field and the 'SAVE BUTTON'.

Drag to Reorder Button and Save

The screenshot displays the GrandCare Designer interface. On the left is a sidebar with navigation options: Accounts, Assessments, Menus, and Logout. The main workspace shows a 'Main Menu' editor. At the top of the editor is an orange header bar with a pencil icon and the text 'Main Menu'. Below this, three menu buttons are displayed on a dark blue background: 'Alex Meal Plan Recipes' (with a calendar icon), 'Photos' (with a photo icon), and 'Alex Meal Plan Recipes' (with a cooking pot icon). A green arrow points to the first 'Alex Meal Plan Recipes' button. Below the menu preview, there are input fields for '* Name or Description' (containing 'Main Menu') and '* Menu ID' (containing 'main'). To the right of these fields is a dropdown for 'Brand' (set to 'Any Brand'). Further right is an orange button labeled 'LAYOUT TYPE: 6 BUTTON'. At the bottom of the interface is a row of action buttons: '< CANCEL', '+ ADD BUTTON', 'EXPORT MENU', 'IMPORT MENU', and 'SAVE MENU'. A green arrow points to the 'SAVE MENU' button.

GrandCare Designer

AUTO-MAVED MENU

Accounts

Assessments

Menus

Logout

Main Menu

Alex Meal Plan Recipes

Photos

Alex Meal Plan Recipes

* Name or Description
Main Menu

LAYOUT TYPE: 6 BUTTON

* Menu ID
main

Brand
Any Brand

< CANCEL

+ ADD BUTTON

EXPORT MENU

IMPORT MENU

SAVE MENU