

How to Add an Assessment Reminder to a Touchscreen

On the Account Communicate page, click on Add Event


The screenshot displays the 'Alex Bedroom' account interface. On the left is a navigation sidebar with options: Home, My Profile, Alerts, Management, Accounts (highlighted), Users, Designer, Enterprise, Groups, Systems, Connectors, Community, Communities, Content, Message Center, Help, and Log Out. The main content area is titled 'Alex Bedroom' and includes a top navigation bar with 'Summary', 'Details', 'Data', 'Caregiving', 'Communicate', and 'Content'. Below this is a toolbar with 'Open Care Menu', 'Show Touchscreen', and 'Video Call'. A prominent orange '+ Add Event' button is located below the toolbar, with a green arrow pointing to it. The main content area also features an 'ALEX BEDROOM CALENDAR' with a date range from Dec 2024 to Feb 2025 and a grid of dates.

Sun	Mon	Tue	Wed	Thu	Fri	Sat
December 2024	December 2024	December 2024	1	2	3	4
5	6	7	8	9	10	11
12	13	14	15	16	17	18

Copyright © 2012-2024 GrandCare Systems LLC
About | Privacy Policy | Help

NOTE: This is a fictional account used for training purposes.

On the Add Event wizard, make sure to select Assessment



The screenshot displays the 'Add Event' wizard interface. On the left is a navigation sidebar with categories: Home, My Profile, Alerts (50), Management, Accounts (selected), Users, Designer, Enterprise, Distributors, Groups, Systems, Connectors, Community, Communities, and Content. The main content area is titled 'ADD EVENT TO ALEX BEDROOM CALENDAR'. It features several form fields: 'Event Type' (a dropdown menu with 'Assessment' selected and highlighted by a green arrow), 'Start Date (required)' (09/12/2024), 'Repeat' (Single - no repeat), 'Start Time' (02:00 PM), and 'Pop-up Tone' (None). A 'Feedback' button is located in the top right corner. The footer contains the text 'Copyright © 2012-2024 GrandCare Systems LLC' and links for 'About | Privacy Policy | Help'.

Complete the remaining fields

ADD EVENT

Feedback

ADD EVENT TO ALEX BEDROOM CALENDAR

Event Type

Assessment

Assessment

Meal Choices

Heart health assessment

COVID-19 assessment

Diabetes assessment

Eating lifestyle assessment

General wellbeing assessment

General lifestyle assessment

Meal Choices

Patient Health Questionnaire (PHQ-9)

Post-seizure assessment

Satisfaction with life assessment

Zones of Regulation

None

Preview Tone

Remind daily until assessment is taken

Home

My Profile

Alerts

Management

Accounts

Users

Designer

Enterprise

Groups

Systems

Connectors

Community

Communities

Content

Message Center

Help

Log Out

grandCARE

Copyright © 2012-2024 GrandCare Systems LLC
About | Privacy Policy | Help


When finished, click Save Event

The screenshot displays the GrandCARE Designer interface. On the left is a navigation sidebar with the following items: Enterprise (highlighted), Groups, Systems, Connectors, Community, Communities, Content, Message Center, Help, and Log Out. A 'Hide' button is located at the bottom of the sidebar. The main content area contains the following fields and controls:

- Start Date (required):** 01/02/2025
- End Date:** Never
- Repeat:** Daily repeat
- Start Time:** 06:00 PM
- Pop-up Tone:** None (with a 'Preview Tone' button below it)
- Remind daily until assessment is taken
- Comments:** (empty text area)
- Buttons:** Cancel and Save Event (highlighted with a green arrow)

At the bottom of the page, the GrandCARE logo is on the left, and the copyright notice 'Copyright © 2012-2024 GrandCare Systems LLC' with links for 'About', 'Privacy Policy', and 'Help' is on the right.

The assessment reminder will display on the designated days and times



Touching the assessment reminder opens the assessment

The screenshot shows a mobile application interface with a dark blue background. On the left side, there is a vertical sidebar with several widgets: a digital clock showing 06:01 pm, a weather widget showing 27°F and Jan 2nd, 2025, three menu buttons labeled 'Main Menu', 'Video Call', and 'Slide Show', and a notification widget showing a bell icon and the text 'No Notifications'. At the bottom of the sidebar are three icons: a gear for settings, a fish for audio, and a question mark for help. The main content area has an orange header bar with the text 'Meal Choices' and a circular icon of a person's head with sound waves. Below the header, the text 'What would you like for breakfast tomorrow?' is displayed. Three large, dark grey buttons with white text are stacked vertically, offering the choices: 'Eggs & Bacon', 'Bagel & Cream Cheese', and 'Pancakes & Syrup'. On the left side of the main area, there is a vertical list of three items: 'Question #1' (highlighted with a yellow arrow), 'Question #2', and 'Question #3'.