


How to Create an Assessment in Designer

From the Account Details page, click on Designer




NOTE: This is a fictional account used for training purposes.

On the Designer Account page, switch to the Assessments tab

The screenshot displays the GrandCare Designer interface. At the top, a dark blue header bar contains a hamburger menu icon and the text "GrandCare Designer". Below this, a light orange sidebar on the left lists four options: "Accounts", "Assessments", "Menus", and "Logout". The "Assessments" option is currently selected. The main content area has a white background. At the top of this area, the text "Alex Bedroom" is displayed. Below it, there are three tabs: "MENUS" (which is underlined in blue), "ASSESSMENTS", and "SETUP". A large green arrow points to the "ASSESSMENTS" tab. The "MENUS" tab is active, showing a list of menu items. At the top right of the menu list, there are two orange buttons: "PUSH MENUS" and "+ ADD". The menu list contains the following items:

Menu Item	Sub-item	Edit Icon	Delete Icon
Games	games - Any Brand		
Health Menu	health - Any Brand		
HEALTHWISE	healthwise - Any Brand		
Main Menu	main - Any Brand		
Message Center	message-center - Any Brand		
Play Menu	play - Any Brand		
Spiritual			

On the Assessments page, click Add



The screenshot displays the GrandCare Designer interface. On the left is a sidebar with navigation links: Accounts, Assessments, Menus, and Logout. The main header area shows the patient name 'Alex Bedroom' and tabs for MENUS, ASSESSMENTS (which is active), and SETUP. Below the header, the title 'Assessments' is centered. To the right of the title are two buttons: 'PUSH ASSESSMENTS' and '+ ADD'. A green arrow points to the '+ ADD' button. Below these buttons is a table listing various assessments.

Assessment Name	Questions	Actions
COVID-19 assessment	COVID-19 Assessment - covid-19 - 7 Questions	✕ ✎
Diabetes assessment	Diabetes Assessment - dia - 5 Questions	✕ ✎
Eating lifestyle assessment	Eating Assessment - eat - 5 Questions	✕ ✎
General lifestyle assessment	Lifestyle Assessment - lif - 6 Questions	✕ ✎
General wellbeing assessment	General Assessment - gen - 3 Questions	✕ ✎
Heart health assessment	Heart Health - chf - 5 Questions	✕ ✎
Patient Health Questionnaire (PHQ-9)	Patient Health Questionnaire - phq9 - 10 Questions	✕ ✎

For the assessment, fill in the fields and add questions

≡

 GrandCare Designer

AUTO-MAVED MENU

 Accounts

 Assessments

 Menus

 Logout

* Assessment Description
Meal Choices

* Assessment ID
meal-choices

* On-Screen Title
Meal Choices

* Assessment Completion Action
Default "Thank You" message

Click *ADD QUESTION* button below to add a question to the assessment

< CANCEL

+ ADD QUESTION

 SAVE ASSESSMENT

For each question, fill in the fields and save

GrandCare Designer

AUTO-SAVED MENU

* Assessment Description
Meal Choices

* Assessment ID
meal-choices

* On-Screen Title
Meal Choices

* Question
What would you like for breakfast tomorrow?

* Question Type
Multiple Choice

Answers

Eggs & Bacon Weight: 1	
Bagel & Cream Cheese Weight: 1	
Pancakes & Syrup Weight: 1	

New Answer Text Weight
1

+ ADD

CANCEL SAVE QUESTION

< CANCEL + ADD QUESTION SAVE ASSESSMENT

When finished, save the assessment

≡

GrandCare Designer

AUTO-MADE MENU

Accounts

Assessments

Menus

Logout

* Assessment Description
Meal Choices

* Assessment ID
meal-choices

* On-Screen Title
Meal Choices

* Assessment Completion Action
Default "Thank You" message

What would you like for breakfast tomorrow?
Multiple Choice - 3 answers

What would you like for lunch tomorrow?
Multiple Choice - 3 answers

What would you like for dinner tomorrow?
Multiple Choice - 3 answers

< CANCEL

+ ADD QUESTION

SAVE ASSESSMENT