


How to Create an Assessment in Designer

From the Account Details page, click on Designer



NOTE: This is a fictional account used for training purposes.

On the Designer Account page, switch to the Assessments tab

The screenshot displays the GrandCare Designer interface. At the top, a dark blue header contains a hamburger menu icon, a user profile icon, and the text "GrandCare Designer". Below the header is a light beige sidebar with four menu items: "Accounts", "Assessments", "Menus", and "Logout". The "Assessments" item is highlighted. The main content area has a white background and features the name "Alex Bedroom" at the top. Below the name are three tabs: "MENU", "ASSESSMENTS", and "SETUP". The "ASSESSMENTS" tab is selected and highlighted with a green arrow. Below the tabs is a large heading "Menus" and two orange buttons: "PUSH MENUS" and "ADD". A table lists several menu items, each with a title, a subtitle, and two icons (edit and delete).

Menu Item	Subtitle	Edit Icon	Delete Icon
Games	games - Any Brand		
Health Menu	health - Any Brand		
HEALTHWISE	healthwise - Any Brand		
Main Menu	main - Any Brand		
Message Center	message-center - Any Brand		
Play Menu	play - Any Brand		
Spiritual			

On the Assessments page, click Add

The screenshot shows the GrandCare Designer interface. The top navigation bar is dark blue with a hamburger menu icon and the text "GrandCare Designer". Below this is a light blue header for "Alex Bedroom" with sub-navigation for "MENUS", "ASSESSMENTS" (which is underlined), and "SETUP". On the left is a vertical sidebar with icons and labels for "Accounts", "Assessments", "Menus", and "Logout". The main content area is titled "Assessments" and contains a table of assessment items. Each item has a title, a subtitle, and two icons (trash and edit). To the right of the table are two orange buttons: "PUSH ASSESSMENTS" and "+ ADD". A green arrow points to the "+ ADD" button.

Assessment Title	Assessment Subtitle	Trash Icon	Edit Icon
COVID-19 assessment	COVID-19 Assessment - covid-19 - 7 Questions		
Diabetes assessment	Diabetes Assessment - dia - 5 Questions		
Eating lifestyle assessment	Eating Assessment - eat - 5 Questions		
General lifestyle assessment	Lifestyle Assessment - lif - 6 Questions		
General wellbeing assessment	General Assessment - gen - 3 Questions		
Heart health assessment	Heart Health - chf - 5 Questions		
Patient Health Questionnaire (PHQ-9)	Patient Health Questionnaire - phq9 - 10 Questions		

For the assessment, fill in the fields and add questions

The screenshot displays the GrandCare Designer interface. On the left is a navigation menu with the following items: Accounts, Assessments, Menus, and Logout. The main area contains a form for configuring an assessment. The form has four fields, each with a green arrow pointing to it:

- * Assessment Description: Zones of Regulation
- * Assessment ID: zones
- * On-Screen Title: Mood Assessment
- * Assessment Completion Action: Default "Thank You" message

Below the form, there is a text instruction: "Click **ADD QUESTION** button below to add a question to the assessment". At the bottom of the screen, there are three buttons: a white "CANCEL" button with a left arrow, an orange "ADD QUESTION" button with a plus sign, and an orange "SAVE ASSESSMENT" button with a save icon. A green arrow points to the "ADD QUESTION" button.

For each question, fill in the fields and save

The screenshot displays the GrandCare Designer interface. On the left is a navigation sidebar with the following items: Accounts, Assessments, Menus, and Logout. The main workspace shows configuration fields for an assessment:

- * Assessment Description: Zones of Regulation
- * Assessment ID: zones
- * On-Screen Title: Mood Assessment
- * Assessment Completion Action: Default "Thank You" message

An "Add Question" dialog box is open in the center, containing the following fields:

- * Question: What zone are you in?
- * Question Type: Zones of Regulation

Three green checkmarks are overlaid on the right side of the dialog box, indicating that the question text, question type, and the "SAVE QUESTION" button have been successfully processed. At the bottom of the main workspace, there are three buttons: a grey "CANCEL" button, a brown "+ ADD QUESTION" button, and a brown "SAVE ASSESSMENT" button with a save icon.

When finished, save the assessment

The screenshot displays the GrandCare Designer interface. On the left is a navigation sidebar with the following items: Accounts, Assessments, Menus, and Logout. The main workspace contains configuration fields for an assessment:

- * Assessment Description: Zones of Regulation
- * Assessment ID: zones
- * On-Screen Title: Mood Assessment
- * Assessment Completion Action: Default "Thank You" message

Below these fields is a preview of the assessment question: "What zone are you in?" with the answer "Zones of Regulation".

At the bottom of the interface are three buttons: CANCEL, ADD QUESTION, and SAVE ASSESSMENT. A green arrow points to the SAVE ASSESSMENT button.